


# POLOS® OIL-LESS VACUUM PUMP


The small-sized POLOS® Oil-less Vacuum Pump is made for continuous run-time, very quiet, and universally applicable. The pump is designed to have no air pollution to the system. Due to its oil-free technology, the POLOS® Oil-less Vacuum Pump can be transported easily and mounted in all positions.




## INCLUDING

- Fitting for tube (6 mm ID, 8 mm OD) - for pump inlet
- 1x Exhaust silencer (with filter) - for pump outlet
- 4x Suction foot mount - for anti-vibration

## DIMENSIONS



## CHARACTERISTICS



A	B	C	A1	B1	D	E	F	H	G1
248	128	186	203	112	186	146	40	M6	1/4-19PT

## TECHNICAL DATA

	Flow	Pressure	Power	Speed	Current	Voltage	Hole	Heat	Sound	Weight
Unit	L/min	Kpa	Kw	Min <sup>-1</sup>	A	V	MM	°C	db(a)	Kg
50Hz	166	-92	0.45	1380	2.21	210/235	6 (ID), 8 (OD)	5-40	58	8.5
60Hz	183	-92	0.52	1450	4.92	110/125	6 (ID), 8 (OD)	5-40	58	8.5